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prismaflex Dose Calculator

Home Filter= M100

Term	Definition
PBP	Pre-Bled Pump
UFR	Ultrafiltration Rate = total replacement + patient weight loss
FF	Filtration Fraction
Effluent	Equals Total replacement + Patient Weight loss + dialysate (if used)
UFR Dosing Delivered	Represents the fluid amounts controlled by the pre, or post, dialyzer, and UFR loss, if they were delivered in a post dialyzer mode. Decreased pre dialyzer clearance effects have been accounted for in all UFR dose calculations to represent the differences of pre/post dialyzer or solute movement, although these effects may have many other variables not present in this tool. Filter size or size, filter pore size, and solute size.

Equations
 $FF = \text{Blood flow based} = \text{Post Replacement} \times 10 / \text{Blood flow}$

Corrected Blood Flow: Prismaflex delivers the Blood flow ordered by increasing the blood pump speed to achieve the ordered blood flow of the access. Prismaflex Corrected Blood flow = Blood Pump speed - PBP*10.

Legend:
 Blue cells: User selected value, or only blue values for different dose results
 Orange cells: Calculated based on values in other cells
 Red cells: Displayed when calculated value exceeded limits

Gambro Renal Products, Inc. has validated the accuracy of the formula for this calculator. However, Gambro is not responsible for arithmetic errors, prescriptions developed in reliance on the user's calculations or patient response to such prescriptions.

Close

GAMBRO
Calculators

Calculate a fluid balance based on patient weight.

% Fluid Accumulation

Calculate a CRRT dose volume based on patient weight.


CRRT Dose Rx

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prismaflex

Leading the way

A graphic of a curved path with four colored spheres (purple, green, pink, and red) on a dark blue background.